

# COMPASS: TEMPEST Discrete DOC Data QAQC

Freshwater Well Test: 2025-11-10

2025-11-18

## Contents

0.1	Run Information . . . . .	1
0.2	Setup . . . . .	2
0.3	Import Data Functions . . . . .	3
0.4	Import Sample Data . . . . .	3
0.5	Assessing standard Curves . . . . .	4
0.6	Assess Check Standards . . . . .	7
0.7	Assess Blanks . . . . .	8
0.8	Assess Duplicates - if there are any . . . . .	9
0.9	Sample Flagging . . . . .	11
0.10	Visualize Data by Plot . . . . .	13
0.11	Convert data from mg/L to uMoles/L . . . . .	15
0.12	Export Processed Data . . . . .	15

## 0.1 Run Information

```
#identify which section you are in  
cat("Run Information")
```

```
## Run Information
```

```
#a link to the Gitbook or whatever protocol you are using for this analysis  
#steph will add this soon
```

```
#anything that needs to be changed do this in the first chunk
```

```
Date_Run = "11/16/25"  
Run_by = "Stephanie J. Wilson"  
Script_run_by = "Stephanie J. Wilson"  
run_notes = " Instrument started, then there was a leak and once that was fixed a new standard curve  
             was run for remaining samples."
```

```
#file path and name for summary file
```

```
raw_file_name = "Raw Data/TMP_AquiferWellTest_NPOC.txt"
#file path and name for the all peaks file
raw_allpeaks_name = "Raw Data/TMP_AquiferWellTest_NPOC_allpeaks.txt"
#file path and name for processed data after QAQC
processed_file_name = "Processed Data/TMP_20251110_FW_WellTest_DOC_Processed.csv"

#check standard concentrations - Update if running different checks:
chk_std_c = 50
chk_std_n = 2

#Log path
Log_path = "Raw Data/COMPASS_TMP_TOCTN_QAQClog_2025.csv"
```

## 0.2 Setup

### 0.3 Import Data Functions

#### 0.4 Import Sample Data

```
## Import Sample Data
```

```
## New names:
```

```
## * ' ' -> '...14'
```

```
## # A tibble: 4 x 4
```

##	sample_name	npoc_raw	tdn_raw	run_datetime
##	<chr>	<dbl>	<dbl>	<chr>
## 1	TMP_AquiferWell_13:30_halffull	0.865	0.165	11/17/2025 10:59:13 AM
## 2	TMP_AquiferWell_13:30_halffull_dup	0.870	0.161	11/17/2025 11:27:18 AM
## 3	TMP_AquiferWell_15:00_full	0.643	0.198	11/17/2025 11:54:54 AM
## 4	TMP_AquiferWell_15:00_full_dup	0.783	0.188	11/17/2025 12:20:54 PM

## 0.5 Assessing standard Curves


```
## Assess the Standard Curve
```

```
## New names:
```

```
## 'geom_smooth()' using formula = 'y ~ x'
```


```
## * '' -> '...18'
```


NPOC Std Curve by Date



```
## 'geom_smooth()' using formula = 'y ~ x'
```

TN Std Curve by Date





```
## [1] "NPOC Curve r2 GOOD"
```

```
## [1] "TN Curve r2 GOOD"
```

## 0.6 Assess Check Standards


```
## Assess the Check Standards
```

```
## New names:
```

```
## * ' ' -> '...14'
```

```
## [1] "Carbon CHECK STANDARD RSD TOO HIGH - REASSESS"
```

```
## [1] "Nitrogen CHECK STANDARD RSD TOO HIGH - REASSESS"
```



```
## [1] ">60% of Carbon Check Standards are within range of the expected concentration"
```

```
## [1] ">60% of Nitrogen Check Standards are within range of the expected concentration"
```

## 0.7 Assess Blanks


```
## Assess Blanks
```

```
## New names:
```

```
## * ' ' -> '...14'
```

```
## [1] ">60% of Carbon Blank concentrations are below the lower 25% quartile of samples"
```

```
## [1] ">60% of Nitrogen Blank concentrations are below the lower 25% quartile of samples"
```



```
## carbon blanks:
```

```
## [1] 0.06378463
```

```
## nitrogen blanks:
```

```
## [1] -0.03622024
```



## 0.8 Assess Duplicates - if there are any

```
## Assess Duplicates
```


```
## # A tibble: 2 x 3
```

```
##   sample_name      npoc_raw_dup tdn_raw_dup
##   <chr>          <dbl>          <dbl>
## 1 TMP_AquiferWell_13:30_halffull 0.870      0.161
## 2 TMP_AquiferWell_15:00_full    0.783      0.188
```

```
##           sample_name npoc_raw tdn_raw      run_datetime
## 1 TMP_AquiferWell_13:30_halffull 0.8650 0.1649 11/17/2025 10:59:13 AM
## 2   TMP_AquiferWell_15:00_full 0.6431 0.1979 11/17/2025 11:54:54 AM
##   npoc_flag tdn_flag npoc_raw_dup tdn_raw_dup
## 1          0.8703      0.1606
## 2          0.7834      0.1876
```

```
##           sample_name npoc_raw tdn_raw      run_datetime
## 1 TMP_AquiferWell_13:30_halffull 0.8650 0.1649 11/17/2025 10:59:13 AM
## 2   TMP_AquiferWell_15:00_full 0.6431 0.1979 11/17/2025 11:54:54 AM
##   npoc_flag tdn_flag npoc_raw_dup tdn_raw_dup npoc_dups_cv npoc_dups_cv_flag
## 1          0.8703      0.1606      0.6465031          YES
## 2          0.7834      0.1876     19.5071025          NO, rerun
##   tdn_dups_cv tdn_dups_cv_flag
## 1    2.776424          YES
## 2    5.562763          YES
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```





```
## [1] "<60% of Carbon Duplicates have a CB <10% - REASSESS"
```

```
## [1] ">60% of Nitrogen Duplicates have a CV <10%"
```

## 0.9 Sample Flagging


## Sample Flagging






## 0.10 Visualize Data by Plot

## Visualize Data





## 0.11 Convert data from mg/L to uMoles/L

## 0.12 Export Processed Data

```
## Export Processed Data
```

```
## # A tibble: 4 x 15
##   Project Experiment Sample_Date Sample_Time Tank_Status Replicate sample_name
##   <chr>    <chr>      <chr>      <chr>      <chr>      <chr>      <chr>
## 1 COMPASS TEMPEST: We~ 2025-11-10 13:30      half full    A          TMP_Aquife~
## 2 COMPASS TEMPEST: We~ 2025-11-10 13:30      half full    B          TMP_Aquife~
## 3 COMPASS TEMPEST: We~ 2025-11-10 15:00      full         A          TMP_Aquife~
## 4 COMPASS TEMPEST: We~ 2025-11-10 15:00      full         B          TMP_Aquife~
## # i 8 more variables: npoc_mgL <dbl>, npoc_uM <dbl>, npoc_flag <chr>,
## #   tdn_mgL <dbl>, tdn_uM <dbl>, tdn_flag <chr>, Analysis_runtime <chr>,
## #   Run_notes <chr>
```

```
#end
```