

COMPASS: TEMPEST Discrete DIC Data QAQC

Freshwater Well Test: 2025-11-10

2025-11-25

Contents

0.1	Import Data Functions	2
0.2	Import Sample Data	2
0.3	Assessing Standard Curves	3
0.4	CRM Check - Don't run chunk if no CRMs run	4
0.5	Assess Check Standards	4
0.6	Assess Blanks	5
0.7	Assess Duplicates	6
0.8	Sample Flagging - Are samples Within the range of the curve?	7
0.9	Visualize Data by Plot	8
0.10	Convert data from mg/L to uMoles/L	9
0.11	Export Processed Data	9

```
##Setup - Change things here & write any notes
```

```
#identify section  
cat("Setup Information")
```

```
## Setup Information
```

```
##### Run information - PLEASE CHANGE  
Date_Run = "11/20/25" #Date that instrument was run  
Run_by = "Stephanie J. Wilson" #Instrument user  
Script_run_by = "Stephanie J. Wilson" #Code user  
run_notes = "Blanks / Background running high, may need to subtract blanks from  
sample concentrations as they were ~3 to 4mg/L. The auto fill IC rxn setting was on  
at the beginning of the run, turned on by a tech during service." #any notes from the run  
samples <- c("TMP") #whatever identifies your samples within the same names  
samples_pattern <- paste(samples, collapse = "|")  
#samples_pattern <- "GCW" #use this instead of the line above if you have only one site code  
chks_name = "Chk_Std_" #what did you name your check standards?  
crm_name = "CRM|crm" #what did you name your CRMS?  
  
##### File Names - PLEASE CHANGE  
#file path and name for raw summary data file  
raw_file_name = "Raw Data/TMP_AquiferWellTest_202511_DIC.txt"  
  
#file path and name for raw all peaks file  
raw_allpeaks_name = "Raw Data/TMP_AquiferWellTest_202511_DIC_allpeaks.txt"  
  
#file path and name of processed data file  
processed_file_name = "Processed Data/TMP_20251110_FW_WellTest_DIC_Processed.csv"  
  
##### Log Files - PLEASE CHECK  
#downloaded metadata csv - downloaded from Google drive as csv for this year  
# Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2025.csv"  
  
#qaqc log file path for this year  
Log_path = "Raw Data/COMPASS_Synoptic_DIC_QAQClog_2025.csv"
```

```
##Set Up Code
```

0.1 Import Data Functions

0.2 Import Sample Data

```
## Import Sample Data
```

```
## New names:  
## * ' -> '...14'
```

```
## # A tibble: 4 x 3  
##   sample_name          ic_raw run_datetime
```

```
##      <chr>                                <dbl> <chr>
## 1 TMP_AquiferWell_1330                    31.4 11/20/2025 9:06:51 PM
## 2 TMP_AquiferWell_1330_dup                 30.8 11/20/2025 9:29:55 PM
## 3 TMP_AquiferWell_1500                    34.1 11/20/2025 9:52:54 PM
## 4 TMP_AquiferWell_1500_dup                 34.5 11/20/2025 10:16:09 PM
```


0.3 Assessing Standard Curves

```
## Assess the Standard Curves
```

```
## New names:
## * ' ' -> '...18'
```


```
## Warning: One or more parsing issues, call 'problems()' on your data frame for details,
## e.g.:
##   dat <- vroom(...)
##   problems(dat)
```

```
## 'geom_smooth()' using formula = 'y ~ x'
```



```
## Warning: Removed 4 rows containing missing values or values outside the scale range
## ('geom_point()').
```

```
## Warning: Removed 4 rows containing missing values or values outside the scale range
## ('geom_line()').
```



```
## [1] "IC Curve r2 GOOD"
```

0.4 CRM Check - Don't run chunk if no CRMs run

```
## Assess the CRMs
```

```
## New names:
## * ' ' -> '...14'
```

```
## [1] "IC crm has a % Difference >25% of expected - REASSESS"
```


```
## Run mean = 29.16
```

```
## Expected = 22.19
```

0.5 Assess Check Standards

```
## Assess the Check Standards
```

```
## New names:
## * ' ' -> '...14'
```




```
## [1] ">60% of IC Check Standards are within range of expected concentration"
```

0.6 Assess Blanks

```
## Assess Blanks
```

```
## New names:
## * ' ' -> '...14'
```

```
## [1] ">60% of Carbon Blank concentrations are lower 25% quartile of samples"
```




```
## carbon blanks:
```

```
## [1] 3.078562
```

0.7 Assess Duplicates


```
## Assess Duplicates
```



```
## [1] ">60% of Carbon Duplicates have a CV <10%"
```


0.8 Sample Flagging - Are samples Within the range of the curve?

```
## Sample Flagging
```



0.9 Visualize Data by Plot

```
## Visualize Data
```

0.10 Convert data from mg/L to uMoles/L

0.11 Export Processed Data

```
## # A tibble: 4 x 12
##   Project Experiment Sample_Date Sample_Time Tank_Status Replicate sample_name
##   <chr>    <chr>      <chr>      <chr>      <chr>      <chr>      <chr>
## 1 COMPASS TEMPEST: We~ 2025-11-10 13:30      half full    A      TMP_Aquife~
## 2 COMPASS TEMPEST: We~ 2025-11-10 13:30      half full    B      TMP_Aquife~
## 3 COMPASS TEMPEST: We~ 2025-11-10 15:00      full         A      TMP_Aquife~
## 4 COMPASS TEMPEST: We~ 2025-11-10 15:00      full         B      TMP_Aquife~
## # i 5 more variables: ic_mgL <dbl>, ic_uM <dbl>, ic_flag <chr>,
## #   Analysis_runtime <chr>, Run_notes <chr>
```

```
#end
```